

ピペットチップ適合表

●=適合

◎=容量に注意

△=未検証

※この適合表をご参照のうえサンプルにてご確認下さい。

		ワトソン																																	
		NEXTYシングル								NEXTYマルチ	NEXTY容量固定				NEXTYオートクレーバブル				NEXTY-S																
容量	品番	NT-2/0.2~2μL	NT-10/1~10μL	NT-20/2~20μL	NT-100/10~100μL	NT-200/20~200μL	NT-1000/100~1000μL	NT-5000/500~5000μL	NTLG-200/20~200μL	NTLG-1000/100~1000μL	NT-810/1~10μL	NT-8200/20~200μL	NT-F10/10μL	NT-F20/20μL	NT-F100/100μL	NT-F200/200μL	NT-F1000/1000μL	NT-F5000/5000μL	NT-A2/0.2~2μL	NT-A10/1~10μL	NT-A20/2~20μL	NT-A100/10~100μL	NT-A200/20~200μL	NT-A1000/100~1000μL	NT-A5000/500~5000μL	NT-S2/0.2~2μL	NT-S10/1~10μL	NT-S20/2~20μL	NT-S100/10~100μL	NT-S200/20~200μL	NT-S1000/100~1000μL	NT-S5000/500~5000μL			
0.1~20μL	204	●	●								●		●							●	●						●	●							
	207	●	●								●		●							●	●						●	●							
	201	●	●								●		●							●	●						●	●							
1~200μL	503			●	●	●			●			●		●	●	●						●	●	●					●	●	●				
	703			●	●	●			●			●		●	●	●						●	●	●					●	●	●				
	705			●	●	●			●			●		●	●	●						●	●	●					●	●	●				
1~250μL	801				●	●			●			●			●	●								●	●					●	●				
10~300μL	505				●	●			●			●			●	●								●	●					●	●				
100~1000μL	502						●			●								●								●									●
	706/804						●			●								●								●									●
	806						●			●								●								●									●
	707						●			●								●								●									●
100~1200μL	805					●			●								●								●									●	
1000~5000μL	401							●										●								●									●
	405																																		
	407																																		
1000~10000μL	409																																		

