

## チップ適合性表

		Sapphire PIPETTE TIPS (STANDARD / LOW RETENTION)							
PIPETTE MODELS	Volume [ $\mu$ l]	10	10 extended	20	100	200	300	1000	1250
	Volume range [ $\mu$ l]	0.2-10	0.5-10	2-20	2-100	5-200	10-300	50-1000	50-1250
	Tip version	Tip/Filter tip	Tip/Filter tip	Filter tip	Filter tip	Tip/Filter tip	Tip/Filter tip	Tip	Tip/Filter tip
	Article no.	7713XX	7723XX	7733XX	7743XX	7753XX	7763XX	7773XX	7783XX
<b>GREINER BIO-ONE</b>									
<b>Single-channel pipettes</b>									
Sapphire	0.2-2 $\mu$ l	●	●						
	1-10 $\mu$ l	●	●						
	1-20 $\mu$ l			●		●			
	10-100 $\mu$ l			■	●	T●/F-			
	20-200 $\mu$ l			■	■	●	●		
	100-1000 $\mu$ l							●	●
	500-5000 $\mu$ l								
	1000-10000 $\mu$ l								
<b>Multi-channel pipettes</b>									
Sapphire	8/12 CH 0.5-10 $\mu$ l	●	●						
	8/12 CH 2-20 $\mu$ l			●		T●/F-			
	8/12 CH 20-200 $\mu$ l			■	■	●	●		
	8/12 CH 20-300 $\mu$ l			■	■	■	●		
<b>BRAND</b>									
<b>Single-channel pipettes</b>									
Transferpette® S	0.5-10 $\mu$ l	●	●						
	2-20 $\mu$ l	■	■	●		T●/F-			
	10-100 $\mu$ l			■		T●/F-			
	20-200 $\mu$ l			■	■	●	●		
	100-1000 $\mu$ l							●	T●/F-
<b>EPPENDORF</b>									
<b>Single-channel pipettes</b>									
Reference®	0.5-10 $\mu$ l	●	●						
	2-20 $\mu$ l			●	●	T●/F-			
	10-100 $\mu$ l			■	●	T●/F-			
	50-200 $\mu$ l				■	●	●		
	100-1000 $\mu$ l							●	●
Research®	0.5-10 $\mu$ l	●	●						
	2-20 $\mu$ l			●	●	T●/F-			
	10-100 $\mu$ l			■	●	T●/F-			
	20-200 $\mu$ l			■	■	●	●		
	100-1000 $\mu$ l							●	●
Research® plus	0.5-10 $\mu$ l	●	●						
	10-100 $\mu$ l			■	●	T●/F-			
	20-200 $\mu$ l				■	●	●		
	100-1000 $\mu$ l							●	●
Research® pro	0.5-10 $\mu$ l	●	●						
	5-100 $\mu$ l				●	T●/F-			
	20-300 $\mu$ l						●		

		Sapphire PIPETTE TIPS (STANDARD / LOW RETENTION)							
PIPETTE MODELS	Volume [µl]	10	10 extended	20	100	200	300	1000	1250
	Volume range [µl]	0.2-10	0.5-10	2-20	2-100	5-200	10-300	50-1000	50-1250
	Tip version	Tip/Filter tip	Tip/Filter tip	Filter tip	Filter tip	Tip/Filter tip	Tip/Filter tip	Tip	Tip/Filter tip
	Article no.	7713XX	7723XX	7733XX	7743XX	7753XX	7763XX	7773XX	7783XX

**Fixed-volume pipettes**

Reference®	50 µl				●	T ● / F -			
	100 µl				●	T ● / F -			
	500 µl							●	●
	1000 µl							●	●
Research®	500 µl							●	●

**Multi-channel pipettes**

Eppendorf	8CH 10-200 µl			■	■	●	●		
Research®	12CH 0.5-10 µl	●	●						
Research® plus	8CH 0.5-10 µl	●	●						
	8CH 10-100 µl			■	●	●	●		
	8CH 20-300 µl			■	■	■	●		

**Electronic pipettes**

Xplorer®	50-1000 µl							●	●
	8CH 15-300 µl			■	■	■	●		

**GILSON**
**Single-channel pipettes**

Pipetman®	P2 (0.5-2 µl)	●	●						
	P10 (1-10 µl)	●	●						
	P20 (2-20 µl)			●		T ● / F -			
	P100 (20-100 µl)			■	●	●			
	P200 (50-200 µl)			■	■	●	●		
	P1000 (200-1000 µl)							●	●
Pipetman® L	P2L (0.2-2 µl)	●	●						
	P10L (1-10 µl)	●	●						
	P20L (2-20 µl)			●		T ● / F -			
	P200L (20-200 µl)				■	●	●		
	P1000L (100-1000 µl)							●	●

**Fixed-volume pipettes**

Pipetman®	F5 (5 µl)			●		T ● / F -			
	F10 (10 µl)			●		T ● / F -			
	F20 (20 µl)			●		T ● / F -			
	F25 (25 µl)				●	T ● / F -			
	F100 (100 µl)						T ● / F -		
	F200 (200 µl)					●	●		
	F250 (250 µl)							●	
	F300 (300 µl)							●	●
	F500 (500 µl)							●	
	F1000 (1000 µl)							●	●

**Multi-channel pipettes**

Pipetman® L	P8X 10L (0.5-10 µl)	●	●						
	P8X 200L (20-200 µl)			■	■	●	●		
	P8X 300L (20-300 µl)			■	■	■	●		

**IKA**
**Single-channel pipettes**

IKA Pette vario	100-1000 µl							●	
-----------------	-------------	--	--	--	--	--	--	---	--

		Sapphire PIPETTE TIPS (STANDARD / LOW RETENTION)							
PIPETTE MODELS	Volume [µl]	10	10 extended	20	100	200	300	1000	1250
	Volume range [µl]	0.2-10	0.5-10	2-20	2-100	5-200	10-300	50-1000	50-1250
	Tip version	Tip/Filter tip	Tip/Filter tip	Filter tip	Filter tip	Tip/Filter tip	Tip/Filter tip	Tip	Tip/Filter tip
	Article no.	7713XX	7723XX	7733XX	7743XX	7753XX	7763XX	7773XX	7783XX
<b>METTLER TOLEDO</b>									
<b>Single-channel pipettes</b>									
Rainin Pipette-Lite™	SL200 (20-200 µl)			■	■	●	●		
Rainin Pipette-Lite™	SL1000 (100-1000 µl)							●	●
<b>NICHIRYO</b>									
<b>Single-channel pipettes</b>									
Nichipet	0.5-10 µl	●	●						
	2-20 µl					T●/F-			
	10-100 µl			■	●	●	●		
	20-200 µl			■	■	●	●		
	100-1000 µl							●	●
Nichipet	8K 30-300 µ				■	■	●		
<b>SARTORIUS</b>									
<b>Single-channel pipettes</b>									
Biohit Proline®	0.5-10 µl	●	●						
Biohit Proline® plus	0.5-10 µl	●	●						
	2-20 µl			●		T●/F-			
	10-100 µl			■	●	●	●		
	20-200 µl			■	■	●	●		
	100-1000 µl							●	●
<b>Multi-channel pipettes</b>									
Biohit Proline®	8CH 0.5-10 µl	●	●						
	8CH 25-250 µl					■	●		
	8CH 50-300 µl				■	■	●		
<b>SOCOREX</b>									
<b>Single-channel pipettes</b>									
Acura® 825	0.5-10 µl	●	●						
	2-20 µl			●		T●/F-			
	5-50 µl				●	T●/F-			
	10-100 µl			■	●	T●/F-			
	20-200 µl			■	■	●	●		
	100-1000 µl							●	●
Calibra® 822	1-10 µl		●						
	2-20 µl			●		T●/F-			
	10-100 µl			■		T●/F-			
	20-200 µl			■	■	●	●		
	100-1000 µl							●	●
<b>Multi-channel pipettes</b>									
Acura® 855	8CH 5-50µl			■	●	T●/F-			
Calibra® 852	8CH 1-10 µl	●	●						
	8CH 20-200 µl			■	■	T●/F-			
	12CH 10-100 µl			■		T●/F-			

		Sapphire PIPETTE TIPS (STANDARD / LOW RETENTION)							
PIPETTE MODELS	Volume [ $\mu$ l]	10	10 extended	20	100	200	300	1000	1250
	Volume range [ $\mu$ l]	0.2-10	0.5-10	2-20	2-100	5-200	10-300	50-1000	50-1250
	Tip version	Tip/Filter tip	Tip/Filter tip	Filter tip	Filter tip	Tip/Filter tip	Tip/Filter tip	Tip	Tip/Filter tip
	Article no.	7713XX	7723XX	7733XX	7743XX	7753XX	7763XX	7773XX	7783XX
<b>STARLAB</b>									
<b>Single-channel pipettes</b>									
ErgoOne®	0.1-2.5 $\mu$ l		●						
	0.5-10 $\mu$ l	●	●						
	10-100 $\mu$ l			■	●	T ● / F -			
	100-1000 $\mu$ l							●	●
<b>Multi-channel pipettes</b>									
ErgoOne®	8CH 30-300 $\mu$ l				■	■	●		
<b>THERMO FISHER SCIENTIFIC</b>									
<b>Single-channel pipettes</b>									
Finnpipette™ F1	1-10 $\mu$ l	●	●						
	2-20 $\mu$ l			●		T ● / F -	●		
	20-200 $\mu$ l			■	■	●	●		
	30-300 $\mu$ l						●		
	100-1000 $\mu$ l							●	●
Finnpipette™ 4500	200-1000 $\mu$ l							●	●
<b>Multi-channel pipettes</b>									
Finnpipette™ F1 MCP8	8CH 1-10 $\mu$ l	●	●						
	8CH 10-100 $\mu$ l			■	●	●	●		
	8CH 30-300 $\mu$ l				■	■	●		
Finnpipette™ 4510	8CH 50-300 $\mu$ l				■	■	●		
Finnpipette™ 4500	8CH 200-1000 $\mu$ l							●	●
<b>VWR</b>									
<b>Single-channel pipettes</b>									
Erg. High Perform.	0.5-10 $\mu$ l	●	●						
	2-20 $\mu$ l			●		T ● / F -			
	10-100 $\mu$ l			■	●	●			
	20-200 $\mu$ l			■	■	●	●		
	100-1000 $\mu$ l							●	●

T = Pipette tip

F = Filter tip

● compatible

- not compatible

■ Pipetting volume limited to the maximum volume of the pipette tips